bodybenefits MASSAGE GUN

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CONAIR
THM001K

bodybenefits MASSAGE GUN

With 4 massage modes, 8 interchangeable heads and 2 hr cordless use, rest, recover & repeat!

READ ALL INSTRUCTIONS BEFORE USING KEEP AWAY FROM WATER

SAFETY PRECAUTIONS

- WARNING:** The polyethylene bags over the product or its packaging may be dangerous. Keep these bags out of the reach of babies and children.
Never use them in cribs, children's beds, pushchairs or playpens. The thin film may cling to the nose and mouth and prevent breathing. A bag is not a toy.
- WARNING:** Do not use this appliance near bathtubs, showers, basins or other vessels containing water.

IMPORTANT:

- Use only the included adapter for recharging the massager.
- Use this massager only for its intended use and as described in this instruction booklet.
- This massager is not intended for commercial use.
- Do not operate under a blanket or pillow. Excessive heating can occur and cause fire, electric shock or injury to persons.
- Do not use any other attachments other than those supplied with this product.
- Never drop the massager or insert any object into an opening.
- Never fall asleep or use while sleeping.
- Do not use outdoors or operate where aerosol (spray) products are being used or where oxygen is being administered.
- Keep the cord away from heated surfaces.
- Do not use the massager if it is damaged. In the event of damage, discontinue use immediately and contact your dealer for instructions on returning it for examination or exchange. If the supply cord is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified person in order to avoid a hazard.

- Never twist or pull on the power cord.
- Always ensure that the voltage to be used corresponds to the voltage marked on the unit.
- For additional protection, the installation of a residual current device(RCD) with a rated operating current not exceeding 30mA is advisable in the electrical current circuit supplying the bathroom. Ask your installer for advice.
- Temperatures sufficiently high to cause burns may occur regardless of the control setting. Do not use on an infant or invalid or on a sleeping or unconscious person. Do not use on insensitive skin or on a person with poor blood circulation. Check the skin in contact with the heated area of the appliance frequently to reduce the risk of blistering.
- This unit is not intended for use by people(including children) with reduce physical, sensory or mental capabilities, or by persons who lack experience or knowledge, unless they are supervised or have been given instructions beforehand on the use of the unit by a person responsible for their safety. Children should be supervised to make sure they do or treat the unit like a toy.

GET TO KNOW YOUR MASSAGER

- Mode 1 : LOW
- Mode 2 : HIGH
- Mode 3 : variable
- Mode 4 : pressure sensor
- On/Off Button with Auto Off : after 10mins for safety and completion of the recommended session time
- Lithium-Ion Rechargeable Batteries 2 hours use time when fully charged
- CORDLESS USE for ultimate autonomy
- Powerful Motor with strong tapping force and 50 percussion per second
- Round Ball Head Attachment for everyday use on larger muscle groups to help broad relief of muscle pain and discomfort

- Bullet Head Attachment use to target deep tissue and tight muscles
- Flat Head Attachment for full body massage
- U-Shaped massage head
- Spiral massage head
- D-Shaped massage head
- Rubber massage head
- Air cushion massage head
- Hard Storage Case stores and protects the percussive massage device and all accessories

CAUTION

- Consult your doctor before using in case of illness or medical condition, including but not limited to pacemaker use, pregnancy, cancer, infection, fracture, or persistent pain.
- Do not use if you have any of the following conditions: diabetes, tuberculosis, benign or malignant tumours, phlebitis or thrombosis, haemorrhages, open or fresh wounds, ulcerated sores, varicose veins, poor circulation or bruised, discoloured, burned, broken, swollen, or inflamed skin or heat sensitivity.
- This product should not be used by any individual suffering from a condition that limits or altogether inhibits the user's capacity to feel or have sensation in any part of the body.
- Do not use this massager on or near your head, neck or throat.
- If you experience pain in a muscle or joint for a prolonged period of time, discontinue use and consult your doctor. Persistent pain could be a symptom of a more serious condition.
- Use of this product should be pleasant and comfortable. Should pain or discomfort result, discontinue use and consult your physician.

NOTE: It is not a medical device for therapeutic purpose .

HOW TO USE

- Before first use, charge the percussive massage device for 5 hours.
- With the device off, install desired attachment head on percussive massage device. Note: Do not change the attachment during the device on.

TIP: For an easy change of the attachment heads, gently pull or insert by turning it from side to side.

- To familiarise yourself with the different speed settings, press the power button on the back of the handle for about 1 second to turn the unit on. The percussive massage device will be in the **Low Speed Setting**. for muscle relaxation, increased blood flow and stress relief. We suggest you start massaging your tight muscles at this slower setting first, without applying much pressure. Each mode is indicated by a blue LED light on the back of the unit. The first blue LED light shows it is in the **Low Speed Setting**.

Note: Do not lean this device on your shoulder or do not place it at an angle on a targeted spot. The accessory heads may be bent or damaged.

- To switch from this low speed setting to the **High Speed Setting:** quick-press the power button again. The massager will now operate at high speed, for a firmer massage, using a stronger tapping force for targeted pain or tension relief. The second blue indicator light on the back of the percussive massage device will light up, indicating it is in **High Speed Setting**.

- If you quick-press the button again, the massager enters the **Variable Speed Setting:** speed varies continuously between low and high speed, mimicking the wave massage technique. The third blue indicator light on the back of the percussive massage device will light up, indicating it is in **Variable Speed Setting**.
- Another quick-press lets you enter the **Pressure Sense Setting:** the harder you press the percussive massager against the targeted area, the more powerful the percussion becomes so you can control the intensity. The fourth blue indicator light on the back of the percussive massage device will light up, indicating it is in **Pressure Sense Setting**.

- To turn off the percussive massage device, press the power button for 1-2 seconds. **The percussive massage device will automatically turn off after 10 minutes of massage time.** If you prefer a longer massage, you can simply turn the unit back on.
- When not in use, your percussive massage device should

be disconnected and stored in a safe, dry location, out of reach of children.

- Use only the included adaptor for recharging the percussive massage device.
- Do not turn on or operate percussive massage device while it is plugged in and being charged.

INDICATOR LIGHTS

- Charger:
 - Blinking blue LED while charging
 - Steady blue LED when product is ready to be used
- Power indication:
 - Blue light indicates there is sufficient power for operation
 - Steady red light indicates that the massager has about 20% of power left
 - Blinking red light indicates that there is less than 10% power left
 - If you see no indicator light, the battery has no power and needs to be charged

NOTE: Do not charge the battery in temperatures below 32 degrees F/0 degrees C.

BODY MASSAGE

Chest:
Massage from centre of chest outward and up toward each shoulder.

Back:
Massage from top of spine to bottom, around shoulder blades and from centre outward.

Abdomen:
Massage from centre outward, in a gentle, circular motion to relax abdominal muscles.

Legs:
Begin at calf, moving up in a circular motion; continue to thighs.

Hands and Arms:
Begin by gently massaging palms of hands and fingers; move upward to forearms in a circular motion, then up to biceps and shoulder muscles.

Feet:
Massage from the heel to the toes in a tight, circular motion.

Do not use this massager on or near the head, neck or throat.

NOTE: This massager is not designed to be used directly on bony areas of the body.

BATTERY DISPOSAL

The THM001K MASSAGE GUN contains two Lithium-ion batteries. To protect the environment, the appliance must be disposed of safely at the end of its useful life. Please take the appliance to a recycling centre, where the internal rechargeable battery should be removed by a professional and recycled separately. For more information about the recycling of electrical and battery operated appliances, please contact your local council office or your household waste disposal service.

ONE YEAR LIMITED WARRANTY

This Conair Appliance is guaranteed for consumer use for one year. Please fill in and return the Warranty Registration as instructed or log into www.conair-asia.com to register. This guarantee is no longer valid in case of alteration or repair by any unauthorized persons and only covers defects under normal use from the date of purchase. If your product does not perform satisfactorily because of defects in materials or manufacture, it will be repaired or replaced under guarantee through the Conair Service Centre. Conair reserves the right to amend terms and specifications without prior notice. To obtain service under this warranty, please send in registration within one month of purchase through either one of the following ways.

THE DATE CODE MENTIONED

WYYXX was mentioned on unit as manufactured date so you can see W as week and Y as year for your reference.

CONAIR KOREA CUSTOMER SERVICE:

Conair Korea Trading Limited
Customer Center: 1899-0209
A/S Center: www.filelec.co.kr

THM001K

IB-20/279

CONTINENTAL CONAIR LIMITED CHINA CREATIVE DEPARTMENT

Date	10/12/2020
Customer	CFEL-KOREA
Model No.	THM001K
Description	IB
Job No.	
Reference No.	IB-20/279
Revision No.	
Color	1C+1C(Black)
Open Size	736mm(W) X 214mm(H)
Folded Size	92mm(W) X 214mm(H)
Material	80gsm Artpaper
Prepared by	Betty bai
Approved by	Ken Xiong / Jay Jeon
Vendor	Chasson
Supplier	

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